# Cavity Liner light cure

## Light-curing radiopaque

### Product description

Cavity Liner light cure is a light-curing, radiopaque one-component material containing calcium hydroxide.

### Fields of application

- Indirect pulp capping.
- As a cavity liner under all restorative materials.

### Application

During application, ensure that the tooth surfaces are perfectly dry and avoid contamination. Instead of using Cavity Liner light cure use, an aqueous calcium hydroxide product on dentine close to the pulp (< 0,5 mm) and for direct pulp capping.

Cavity Liner light cure can be applied directly from the syringe.

Cavity linings with a thickness in excess of 1 mm need to be built up in layers and each layer cured using blue light for at least 20s.

If increased adhesion values are required (e. g. when using Cavity Liner light cure for adhesive restorative treatments), a suitable adhesive can be used prior to the application of Cavity Liner light cure. To this end, follow the adhesive manufacturer’s instructions for use.

Carefully remove any excess material prior to light-curing using suitable tools. Do not remove polymerised excess material with a probe but, instead, carefully use a rotating instrument and then place the restoration according to the manufacturer’s instructions.

### Information, precautionary measures

Do not use Cavity Liner light cure in the restoration margins. In the marginal area, the permanent restoration should always be placed on caries-free dental hard tissue which, when using the adhesive technique, must have been conditioned accordingly.

Cavity Liner light cure adheres to methacrylate-based products, e. g. crown and bridge materials, temporary sealing materials and certain impression materials (e. g. polyether), which can lead to removal of the cavity lining. In such cases, application of a separating agent is recommended (e. g. Vaseline, glycerine gel etc.). Remove the insulating material prior to placement of the final restoration.

Once dispensed, Cavity Liner light cure must not be exposed to direct light (operating light, sunlight) and must be applied as quickly as possible.

Cavity Liner light cure contains methacrylates, BHT and amines. Cavity Liner light cure should not be used in case of known hypersensitivities (allergies) to any of these ingredients. Calcium hydroxide has a caustic effect on the skin and mucosal tissue, it must not be allowed to come into contact with the eyes!

Phenolic substances, especially products containing eugenol and thymol, impede the curing of Cavity Liner light cure. The use of zinc oxide eugenol cements or other eugenol-containing materials in combination with Cavity Liner light cure should therefore be avoided.

Hold the output window of the light as close to the surface of the restoration as possible, as curing may otherwise be insufficient. Inadequate curing may lead to loss of strength and to discomfort.

### Storage

Store at temperatures of 4° C to 23° C protected from light and moisture.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and / or written information has been given to the best of our knowledge and without obligation. Our information and / or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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MP Medical Product GmbH