DIRECTIONS FOR USE**GB****17% EDTA solution**

Chelating agent for removing the smear layer and opening dentinal tubules in root canal irrigation.

Composition

Tetrasodium edetate (EDTA) 17%, cetrimonium bromide 0.9%, purified water ad 100%.

Properties

17% EDTA solution helps soften the dentine walls of the root canal and remove the smear layer in root-canal irrigation. Its effective mechanism is the chelation of calcium ions and the removal of inorganic debris.

Indication

Irrigation, removal of debris and rinsing of the root canal.

Directions for use

Irrigate the canal with the solution using e.g. a syringe and allow to soak for 1 minute before rinsing thoroughly with 3% sodium hypochlorite solution. Aspirate excess liquid, then dry the canal.

Warning

Avoid extrusion of the solution beyond the root apex. When used alternately with sodium hypochlorite, thoroughly rinse the root canal, EDTA reduces the concentration of active chlorine from the hypochlorite solution.

In case of skin contact rinse thoroughly with water and soap. If skin irritation occurs, consult a physician.

In case of contact with eyes, rinse immediately with plenty of water for 15 minutes while holding the eyelid open and seek medical advice immediately.

If swallowed, rinse mouth with water and seek medical advice immediately. Do not use in known allergies to any of the product ingredients.

Storage

Reseal the vial firmly immediately after use. Store at room temperature (between 15°C and 25°C). Protect from light. If the product is stored at temperatures below 15°C, small EDTA crystals may form. These crystals will disappear once the temperature goes above 15°C again.

Packaging

1 vial (15 ml)

REF

11565

For dental professional use. Made in Switzerland.

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