
FIBRAPOST LUX

Light conducting Glass Fibre reinforced Composite Posts for the coronal-radicular dental reconstruction

Directions for use

This instruction manual refers to the use of SEALACORE DC (PD; Switzerland) products. **Note:** for the use of another bonding and/or core build up system read the instructions of the manufacturer.

Endodontically prepare the canal in the usual manner, clean and dry. Obturate the apical portion.

- 1 Choose the post size and measure the length of the canal.
- 2-3 Drill the canal with the corresponding Drill-A-Post (800-1'200 RPM), placing its depth stopper at the required length.
- 4 Check the size of the post and its position in the canal.
- 5 Excess in length may be adjusted with a diamond disc. **Carry it out outside of the mouth**, under aspiration. Clean the post with alcohol.
- 6 SEALACORE DC Nano-reinforced Self-Etch Bond. Mix 1 drop of Liquid A and 1 drop of Liquid B on a mixing pad with the application brush for approximately 2 seconds. Avoid exposure to intensive light. Apply a layer with the application brush to the enamel/dentine and rub into the tooth surface for 20 seconds.
- 7 Dry the adhesive layer for at least 5 seconds with an air syringe.
- 8 Do not light cure at this step of the procedure.
- 9 Apply SEALACORE DC Composite directly into the root canal using thin intraoral tips. **Note: Do not light cure at that stage!**
- 10 Insert the FIBRAPOST making sure to leave a slight coronal excess for further anchorage of the core.
- 11 Light-cure for at least 40 seconds.
- 12 **Core build up.** Place a matrix around the prepared tooth and apply SEALACORE DC Composite directly from the cartridge/double-barrel syringe. Place the tip of the polymerization lamp as close as possible to the surface of the composite and polymerize with halogen light for 40 seconds. In case of a larger core build-up, generally apply in layers of 2 mm and cure each layer for 10 seconds.
Final Polymerisation: Light cure 40 seconds or wait 5 min. After curing, the composite can be worked immediately.

Properties

FIBRAPOST LUX are made from specially treated Glass Fibres (approx. 65 %) reinforced with, and firmly bound to, an UDMA (Urethane Dimethacrylate) resin (approx. 20%) known for its resistance, safety and capacity to bind to Composite resins used for post cementation.

Features and benefits

- Light transmission, both axially (towards the apex) and laterally, by diffusion
- Optimal elastic modulus (~ 35 GPa) close to the one of the dentine
- High resistance (~ 1400 MPa) thanks to interlaminar fibres cohesion
- Translucent shade for aesthetic reconstruction
- Convenient Radioopacity (~ 300% Al (Ø 2.5 mm))
- High adhesion capacity between post and cement
- Biocompatibility

Indication

Aesthetic post-retained restorations and core building-up of endodontically treated non-vital teeth.

Sizing and shape

4 sizes according to canal size allowing the preservation of maximum healthy dentine. Cylindro-conical shape.

Contraindications

- Do not use Fibrapost:
- In case of allergy to methacrylates
- If the root is not well anchored in bone
- In the mesial root of a mandibular molar
- In the buccal root of a maxillary molar
- In case of insufficient tooth structure
- In teeth with ferrule less than 2 mm
- Fibrapost must not be sterilized before use

Disinfection

Disinfect with alcohol before using. Single use.

Storage Conditions

5°C – 30°C. Protected from humidity.

Product presentation in blister

Do not add any liquid or other material into the blister.

Packaging

Assortment classic
6 posts each N° 1 to 4 = 24 posts
1 reamer each N° 1 to 4 = 4 reamers

REF
22589

Refills (in racks of 12 pieces)

N° 1	22590
N° 2	22591
N° 3	22592
N° 4	22593
N° 1 to 4 (3 pieces of each size)	22594

Refills in blister (box of 5 pieces)

N° 1	22641
N° 2	22642
N° 3	22643
N° 4	22644

Size	N° 1	N° 2	N° 3	N° 4
(L) Total length in mm	19	19	19	19
(l) Length of the tapered part in mm	11	11	11	11
(D1) Ø in mm at the tip	0.85	1.06	1.26	1.47
(D2) Ø in mm at 11 mm	1.30	1.50	1.70	1.90

For dental professional use. Made in Switzerland.

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