DIRECTIONS FOR USE**Irriflex®**

Endodontic needle for root-canal irrigation.

Indication

Root canal irrigation.

Properties

Flexible cannula with radial opening for liquid ejection along the canal walls.

Direction for use

1. Clean the packaging by pulling along pre-cut notch provided for this purpose.

2. Prepare the syringe with one of the compatible irrigation solutions: stabilized sodium hypochlorite up to 6%, EDTA 17%, Chlorhexidine digluconate 2%, Acid citric up to 40%. Dual Rinse® HEDP. Make sure the irrigation solution has been properly stored to avoid the formation of any crystal which could block the tip of the needle.

3. Insert the "Secure lock + part of the needle" into the male nose of the previously filled Luer-Lock syringe and lock the needle until it clicks into place.

4. Make sure the needle is firmly locked on the syringe and expel any residual air from the syringe.

5. Identify the 4 depth marks on the canal (18mm, 19mm, 20mm and 22mm). A stopper may be used to define the working length.

6. Before entering the canal, check the pressure needed on the syringe barrel to create two equal jets exiting the needle. This step also allows to control the flowability of the irrigation solution.

Irrigation

7. Start irrigating in the pulp chamber and gently insert the cannula into the canal orifice. Pre-bend the tip to facilitate access into posterior teeth.

8. Use a stopper to indicate a distance of 1 mm from the canal orifice.

9. Determine the working length with 1 or 2 ml.

10. Use your shaping sequence. As a rule, after each active instrument used in the sequence, irrigate with 1 or 2 ml of solution, check patency and re-irrigate to eliminate debris. Make sure that the cannula can move freely to ensure constant flowability of the irrigation solution.

11. As a final rinse, it is recommended to flush the canal with EDTA followed by sodium hypochlorite.

12. After use, disconnect the cannula and dispose according to local regulations.

13. Check the length of the canal to be irrigated and use the 18-, 19-, 20- and 22-mm length markings to ensure proper positioning of the needle. Make sure that the needle can move freely to ensure constant flowability of the irrigation liquid.

Compatible Irrigation solutions

Soda hypochlorite up to 6%, EDTA 17%, Chlorhexidine digluconate 2%, Acid citric up to 40%, Dual Rinse® HEDP.

Adverse reactions

None known.

Contraindications

None known.

Warmer messages

Single-use device. Do not reuse.

Storage

Keep away from direct light.

Shelf-life

5 years. Do not use the product after the expiry date.

Disposal

Discard each cannula as per local regulation.

Désinfection, cleaning and sterilisation

This product has been manufactured under clean-room conditions. It is ready for use. Do not sterilise.

Product description

Box of 20 needles.

Box of 40 needles.

Needles are individually packaged in plastic pouch.

For dental professional use. Made in Switzerland.

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