DIRECTIONS FOR USE

Zinc Oxide

Co-formative agent of ZOE cement

Composition

Zinc Oxide 100%

Properties

When mixed with Eugenol, Zinc Oxide powder forms a cement (paste) used in dentistry.

Indications

- ▶ Temporary obturation
- ▶ Indirect pulp capping

Direction for use

Mix the zinc oxide powder with eugenol in the proportions of powder / liquid (P/L) between 1.1 and 1.7 (g/g).

Setting starts in the presence of moisture. An increase in humidity or temperature accelerates the setting. The setting time (tested according to ISO 3107) is 6h 50-16h according to the P/L ratio used.

EN

Warning

In case of contact with the skin wash with plenty of soap and water. In case of contact with the eye rinse cautiously with soft, clean water for 15 minutes, holding the eyelids open, and call a doctor immediately. In case of swallowing rinse the mouth with water, drink a lot of water and consult a doctor. Keep the person exposed at rest, do not force vomiting. Do not use in cases of known hypersensitivities to any of these ingredients.

Storage

Store at room temperature between 5°C and 30°C.

Packaging

	REF
Bottle of 100g	11431
Bottle of 250g	11432
Bottle of 500g	11430
Bottle of 1000g	11433

For dental professional use.

PD 200617

Zinc Oxide

Powder . Regular setting

Blended into a paste with Eugenol, Zinc Oxide powder becomes a temporary cement which is also used for indirect pulp capping.

For dental professional use . Made in Switzerland

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