

















	EN Symbols	FR Symboles	DE Symbole	IT Simboli	ES Símbolos	PT Símbolos
	Manufacturer	Fabricant	Hersteller	Produttore	Fabricante	Fabricante
	European Authorised Représentative	Mandataire autorisé dans l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft	Rappresentante autorizzato nella Comunità Europea	Representante autorizado europeo	Representante Europeu Autorizado
	Authorised representative for Switzerland	Mandataire établi en Suisse	Bevollmächtigter in der Schweiz	Mandatario in Svizzera	Representante autorizado para Suiza	Representante autorizado para a Suíça
	Date of manufacture	Date de fabrication	Datum der Herstellung	Data di produzione	Fecha de fabricación	Data de fabrico
	Use by date	Date limite d'utilisation	Verwendbar bis	Utilizzare entro la data indicata	Fecha de caducidad	Utilizar até
	Batch number	Numéro de lot	Chargennummer	Codice di lotto	Número de lote	Número de lote
	Catalogue number	Référence catalogue	Katalognummer	Numero di catalogo	Número de referencia	Número de referência
	Non sterile	Non-stérile	Nicht steril	Non sterile	No estéril	Não estéril
	Sterilizable in a steam sterilizer (autoclave) at temperature specified	Stérilisable à la vapeur (autoclave) à la température spécifiée	Sterilisierbar in einem Dampfsterilisator (Autoklav) bei angegebener Temperatur	Sterilizzabile in uno sterilizzatore a vapore (autoclave) alla temperatura specificata	Esterilizable en un esterilizador de vapor (autoclave) a la temperatura especificada	Esterilizável num esterilizador a vapor (autoclave) à temperatura especificada
	Fragile; handle with care	Fragile ; manipuler avec soin	Zerbrechlich; mit Sorgfalt handhaben	Fragile; maneggiare con cura	Frágil; tratar con cuidado	Frágil; manusear com cuidado
	Keep away from sunlight	Conserver à l'abri de la lumière du soleil	Von Sonnenlicht schützen	Conservare al riparo dalla luce solare	Mantener alejado de la luz del sol	Manter afastado da luz solar
	Keep away from rain	Craint l'humidité	Trocken aufbewahren	Conservare a secco	Mantener alejado de la lluvia	Manter afastado da chuva
	Temperature limit	Limite de température	Temperaturbegrenzu ng	Limite di temperatura	Límite de temperatura	Limite de temperatura
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	Non riutilizzare	No reutilizar	Não reutilizar
	Consult instructions for use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso	Consultar instrucciones de uso	Consulte as instruções de utilização
	Caution	Attention	Vorsicht	Attenzione	Precaución	Cuidado



	Green dot	Point vert	Grüner Punkt	Punto Verde	Punto verde	Ponto verde
	CE marking of conformity	Marquage de conformité CE	CE-Kennzeichnung der Konformität	Marchio di conformità CE	Marcado CE de conformidad	Marcação CE de conformidade
	Medical device	Dispositif médical	Medizinprodukt	Dispositivo medico	Producto sanitario	Dispositivo médico
	Unique Device Identifier	Identifiant unique des dispositifs	Einmalige Produktkennung	Identificativo unico del dispositivo	Identificador único del producto	Identificação única do dispositivo