

The original Tofflemire

Matrix Retainers . Porte-matrices . Matrizenspanner
Portamatrici . Porta-matrices . Porta matriz



Stainless Steel

Junior
+12 Matrix Asstd.

REF

41504

41501

= (41504 + 41528)



Senior
+12 Matrix Asstd.

41505

41502

= (41505 + 41528)



Universal
+12 Matrix Asstd.

41506

41503

= (41506 + 41528)



• VASELINE . VASELINA . VASELINA

• FOR SMOOTH WORKING . POUR UNE PARFAITE MANIABILITE
• FÜR EINE ANGENEHME BENÜTZUNG
• PER MIGLIORARE L'IMPIEGO . PARA UNA UTILIZACIÓN SUELTA
• PARA UMA MELHOR UTILIZAÇÃO

Various possibilities of use . Nombreuses possibilités d'emploi
Vermehrte Anwendungs möglichkeiten . Svariate possibilità di applicazione
Varias modas de empleo . Muitas possibilidades de utilização
Disinfection and sterilization with each use . Désinfection et stérilisation
avant chaque usage . Desinfektion und Sterilisation mit jedem Gebrauch
Disinfezione e sterilizzazione con ogni uso . Desinfección y
esterilización . Desinfecção e esterilização com cada uso con cada uso




Tofflemire Matrix Bands

Matrices . Matrizen . Matrici . Matrices . Matrizes

Alloy Stainless Steel

Thickness Thin 0.05 mm, x-thin 0.035 mm

Actual sizes



Pack of 12 pieces	REF	REF
Adults	thin	x-thin
#1	41521	41511
#2	41522	41512
#3	41523	41513
#11	41524	41514
Children		
#13	41525	41515
#14	41526	41516
#15	41527	41517
Assorted		
4× #1		
3× #2	41528	41518
1× #3-15		

 **Produits Dentaires SA**
Rue des Bosquets 18
1800 Vevey . Switzerland



Made in Switzerland
pd-dental.com

For dental professional use . Pour usage dentaire professionnel
Für den professionellen zahnärztlichen Gebrauch . Per uso dentale professionale
Para uso dental profesional . Para uso dental profesional